LC-Base plate

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Instructions for use (EN)

Read and keep before using the product! Observe the safety instructions!

General safety instructions

Observe safety data sheet!

H317: May cause an allergic skin reaction. Avoid contact of uncured material with skin/mucous membrane and eyes. It may be slightly irritating and in rare cases may cause sensitization to Methacrylates.

Use protective equipment: Safety glasses, suitable protective gloves with protection function to Methacrylates (e.g. Polyethylene-gloves, 4h-gloves). Do not wear contaminated work clothing outside the workplace.

After skin contact: Wash immediately with plenty of soap and water. Medical treatment necessary if symptoms occur (irritation of the skin).

After eye contact: Immediately flush with running water for 10 to 15 minutes with the eyelid open and consult an ophthalmologist.

When grinding the polymerized sheets, use a suitable dust extractor and follow the instructions for use of the abrasive manufacturer, especially with regard to dust protection. If necessary, use a self-contained breathing apparatus. All serious incidents related to the product must be reported to the trademark owner and the competent authority of the member state where the user and/or patient is established. Please direct any feedback regarding the product to: registration@mueller-omicron.de

General description: LC-base plate is a pre-shaped, modelable, light-curing composite for the production of individual base plates in the dental laboratory or dental practice. It features low shrinkage, low water absorption and high dimensional stability. The thickness of the plate is about 1.4 mm. The product is available in the color pink.

Composition: Mixture of different acrylates / methacrylates, fumed silica, glass beads, photochemical initiators, vaseline.

Intended use

Warnings: LC-base plate is light-sensitive, so that polymerization can already be triggered by exposure to daylight. Therefore, the material should not be exposed to light unnecessarily (removal from the can only immediately before adaptation to the model, immediate sealing of the can after removal).

Intended use:

LC-base plate are for dental use only: for production of a occlusal base plate / for intraoral support pin registration. The user is obliged to examine the material for suitability and usage possibilities for the intended purposes on his own responsibility before using it, especially if those are not listed in this instruction manual. The manufacturer accepts no liability for damage resulting from any other use.

Clinical benefit: The clinical benefit of light-curing base plates lies indirectly in the replacement of the dental apparatus and thus in a restoration of masticatory function, since precision and/or functional imaging are essential for the accurately fitting fabrication of dentures.

Indications: Bite registrations and Support pin registrations

Contraindication: In case of known hypersensitivity to methacrylates, skin contact should be avoided.

User profile: LC-base plate is a product for professional use by dental technicians or dentist. It is intended for use in the dental laboratory or dental practice. No special training is required. The finished products made of the material (base plates, support pin registrations) are intended for professional use by dentists in a dental practice. Special training is not required.

Patient group: The use of LC-base plate is not limited to any particular patient group. Contraindications must be observed.

Undesirable side effects: LC-base plate contains methacrylates. These may cause allergic skin reactions. In rare cases sensitization is possible.

Interactions: No interactions are known.

Special storage instructions: Keep the can tightly closed. Store protected from light and at 4 - 25°C.

Shelf life: 3 years from manufacture. Do not use after the expiry date stated on the label.

Waste treatment: Do not allow product or primary packaging which has not been emptied of all residues to reach sewage system, water, subsoil or soil. Disposal according to official regulations. Non-contaminated and completely empty packaging as well as fully cured products can be disposed of with normal municipal waste (household waste). Observe waste separation.

Application

- The situation model should be prepared as required and insulated with an alginate insulation. Undercutting areas of the model should be blocked out before.
- Open the can at the handle tab, remove one LC-base plate incl. protective fleece and immediately reseal the can tightly.
- The protective fleece is removed from the LC-base plate sheet.
- The sheet is adapted to the insulated model and trimmed at the edges.
- The sheet is then polymerized on the model in a blue light or UVA light device for 5 minutes. The instructions for use of the light unit used must be observed!
- The sheet is removed from the model and polymerized from the back for another 5 minutes.
- The individual base plate made in this way can be processed immediately.
- The inhibition layer (adhesive layer on the top of the material) allows subsequent correction (e. g. reinforcement of the sheet). The inhibition layer can be easily removed with a moistened tissue with an alcoholic cleaning agent.

Packing form: REF 0000LC018, box with 50 plates, colour pink

Symbols used in product labeling:

- | | | | | | | | | | |
|--|---|--|-------------------------------------|--|------------------------------|--|-------|--|--------------|
| | Manufacturer | | Medical device Item | | number | | Batch | | Usable until |
| | Protect from sunlight | | Temperature limit | | Observe instructions for use | | | | |
| | Attention! Observe safety / warning notes, precautionary measures in instructions for use | | May cause an allergic skin reaction | | | | | | |

A Quality Brand from

Müller-Omicron GmbH & Co. KG,

Schlosserstr. 1, 51789 Lindlar / Germany, phone +49(0)2266-47420, www.mueller-omicron.de

MEGADENTA Dentalprodukte GmbH Carl-Eschebach-Str. 1A, D-01454 Radeberg, Germany

REV 00-2021