



Dentodril® ID liquid

**Universal disinfectant for
rotating instruments
READY-TO-USE**

Advantages:

- ▶ Broad spectrum of activity
- ▶ Very good material compatibility
- ▶ Longer lifetime of the drills due to mild pH-value
- ▶ With corrosion protection

Product Type:

Highly effective and cleaning active, ready-to-use combination preparation with comprehensive spectrum of activity against bacteria (incl. TBC), yeasts and enveloped viruses. The preparation contains an effective corrosion protection. VAH-listed.

Type of application:

For cleaning and prophylactic disinfection of dental drills and rotating instruments, such as carbide drills, diamond cutters, ceramic abrasives, rubber polishers, root canal instruments and other stainless steel instruments in an ultrasonic bath or bur bath.

Indication for use:

Always use Dentodril® ID liquid undiluted! Immediately after use, place the instruments in the ready-to-use solution; all parts to be disinfected must be completely covered. Remove all parts from the bath at the earliest after the exposure time has expired, rinse with water and dry. Do not mix with other disinfectants or cleaners, not suitable for aluminium. It is recommended to wear protective gloves and eye protection.

Spectrum of effectiveness:

Bactericidal (incl. MRSA) as per EN 13727, EN 14561, VAH; levurocidal (Candida Albicans) as per EN 13624, EN 14562, VAH; Mycobactericidal (M. terrae and M. avium) as per EN 14348, EN 14563; "limited virucidal" (enveloped viruses such as SARS-CoV-2, BV, HIV, HCV) as per DVV/RKI and EN 17111.

Application times:

Bactericidal (incl. MRSA) as per EN 13727, EN 14561, VAH (dirty conditions): conc.: 5 min.

Levurocidal (Candida Albicans) EN 13624, EN 14562, VAH, (dirty conditions): conc.: 5 min.

Mycobacteria (M. terrae and M. avium) as per EN 14348, EN 14563 (clean conditions): conc.: 5 min.

"Limited virucidal" (enveloped viruses e.g. SARS-CoV-2, HBV, HIV, HCV) as per DVV/RKI and EN 17111 (dirty conditions): conc.: 5 min.

Ingredients (100 g contains):

0,396 g didecyldimethylammonium chloride, 0,2475 g N-(3-aminopropyl)-N-dodecyl-propane-1,3-diamine.

Chemical-physical data:

- ▶ Appearance: colourless
- ▶ Density: 1.01 g / cm³
- ▶ pH value: 9.8

Safety notes:

H 412 Harmful to aquatic organisms, with long-term effects. Avoid release into the environment.

Notes on disposal:

Disposal of the contents/container in accordance with local/regional/national/international regulations. Detergents VO 648/2004/EG: Contains < 5 % non-ionic surfactants, < 5 % cationic surfactants, fragrances.

Delivery forms:

Packing unit	REF
5 litre canister	73005000

Particularly suitable accessory:

Product:	REF
Canister cock	701000911



Shelf life:

3 ½ Years from the date of manufacture.

Contact:

Further information about reports / test results and about the 360° hygiene program is available on request.

Tel.: +49 (0) 2266 4742 - 0

E-Mail: info@mueller-omicron.de

A Quality Brand from

Müller-Omicron GmbH & Co. KG

Schlosserstr. 1

D-51789 Lindlar

Tel.: +49 (0) 2266 4742-0

Fax: +49 (0) 2266 4742-23

www.mueller-omicron.de

www.mueller-omicron.com