



Disinfectant for Impression CONCENTRATE

Your advantages:

- ▶ Broad spectrum of activity including TB
- ▶ Contains no aldehydes, alcohols or toxic ingredients
- ▶ Does not change dimensions

Preparation type:

Concentrate for disinfection of dental impressions.

Fields of application:

For disinfection of all dental impression materials.

Indication for use:

Dilute the concentrate with 10°C to 30°C warm water (at least drinking water quality) (use dispensing bottle). Place the respective impression in the ready-to-use disinfectant solution; the impression must be completely covered with the solution. Especially in the case of alginates and hydrocolloids, the application times must be observed (remove them immediately once the application time is over). After removal, the impression should be rinsed with water. Do not mix with cleaning agents or other disinfectants! The working solution can be disposed of in the sewage water system. The dosage quantities and application times should be observed as precisely as possible and should not be below these values. For professional users only.

Spectrum of effectiveness + application times:

Bacteria as per EN 13727, EN 14561,
VAH (Association of applied Hygiene)
(dirty conditions): 2.0 % 15 min.

Candida albicans as per EN 13624,
EN 14562, VAH (dirty conditions): 2.0 % 15 min.

„Limited virucidal effect“ (enveloped viruses such as
SARS-CoV-2, HBV, HIV, HCV) as per DVV/RKI and EN
17111 (dirty conditions): 2.0 % 15 min.

TB (M. terrae) as per EN 14348
and EN 14563 (clean conditions): 4.0 % 30 min.
3.0 % 60 min

Ingredients:

2.75g N-(3-aminopropyl)-N- dodecylpropane-1.3-
diamine;
4.40g N,N-didecyl-N,N-dimethylammonium chloride



Chemical-physical data:

- ▶ Appearance: green
- ▶ Density: 1,01 g / cm³
- ▶ pH value at 20° C: > 10

Special notes / safety instructions:

Safety instructions: H 302 Harmful if swallowed. H 314 Causes severe skin burns and eye damage. H 410 Very toxic to aquatic life with long lasting effects. Avoid release to the environment. Wear protective gloves/eye protection. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower]. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Immediately call a POISON CENTER/doctor.

Notes on disposal:

Dispose of the contents/container in accordance with local/regional/national/ international regulations. Contains didecylmethylammonium chloride, N-(3-aminopropyl)- N-dodecylpropane-1.3-diamin.

Delivery forms:

Packing Unit	REF
2 litre bottle	74102000

Particularly suitable accessory:

Product	REF
1 L Universal metered-dose-squeeze bottle, empty	701001000



Shelf life: 3½ Years from the date of manufacture.

Contact:

Further information about reports / test results and about the 360° hygiene program is available on request. Tel.: +49 (0) 2266 4742 - 0
E-Mail: info@mueller-omicron.de

A Quality Brand from

Müller-Omicron GmbH & Co. KG
Schlosserstr. 1
D-51789 Lindlar
Tel.: +49 (0) 2266 4742-0
Fax: +49 (0) 2266 4742-23
www.mueller-omicron.de
www.mueller-omicron.com

